

couture

recliner



 **LINET**[®]
AMERICAS

couture recliner

Thoughtfully designed and engineered for use in rigorous 24/7 healthcare environments, the Couture recliner from LINET certainly lives up to it's name.

Generously padded all over, and with a gas cylinder mechanism allowing infinite recline positioning, the Couture recliner is a treat to relax in.


Equally important to patient comfort is ease of operation for the caregiver. When developing the Couture Recliner we gave extra consideration to this, which is evidenced in the features below.

Optional arm cap is available in a choice of wood or polyurethane.



* Shown with optional wood arm cap, central locking system and Trendelenberg feature. Cover image also features optional IV pole, side tray and foley bag holder.

features



Optional "Central Locking System" allows all four casters to lock simultaneously



Optional flip-up side tray can be specified to be interchangeable from left to right side of the recliner



Optional slide out foot rest provides easier patient transport



Optional Trendelenberg feature is easily activated from either side of the recliner




Ergonomically placed handles enable easier patient transport


Wood arm finish selection:




Hardrock Maple (HM)




Candlelight (CL)




Rustic Cherry (RC)




Summer Flame (SF)



Mahogany (MA)




Cherry Blossom (CB)



Chocolate (CH)

Now Available!

The Couture Recliner now features an optional "swing-up" arm for easier patient transfers. Innovative design features a pneumatic spring and is efficient, safe and easy to use.



LINET Americas, Inc., 10420-R Harris Oaks Blvd, Charlotte, NC 28269
Phone: 877.815.9895 | Fax: 704.248.5655 | www.linetamericas.com

