

# Arezzo 3-Position Sleeper

The innovative Arezzo Sleeper provides an attractive 3-position seating/sleeping solution for patient rooms. It can be specified to open to the left or the right, so it's ideal for small spaces unable to accommodate typical convertible furniture that pulls out into the room.

## Healthcare Features

The Arezzo Sleeper incorporates the following features specific to healthcare seating products:

- Sleep surface is 28" x 80"
- Heavy-duty all steel frame construction
- Central locking caster option engaged with foot pedals
- Field-replaceable seat and bolster foam
- Field-replaceable seat, bolster and end frame upholstery
- Polyurethane arm cap option
- Replaceable wood arm cap option
- Seat and back foam meets or exceeds CAL 117-2013 flammability requirements
- Optional CAL 133 compliance
- Optional Moisture Barrier
- Limited lifetime warranty



## Construction Specifications

Steel end frame assemblies are constructed of powder coated 14-gauge 1-5/8" tubular steel. The wheels are made of long-lasting neoprene.

The sleeper can be specified with an upholstered arm or with European steamed beech or black polyurethane arm caps. It can be fully upholstered or with wood front panel and optional storage drawer. The arm caps and front panel are available in our 7 standard wood finishes:

- |                        |                       |                  |
|------------------------|-----------------------|------------------|
| • Candlelight (CL)     | • Cherry Blossom (CB) | • Chocolate (CH) |
| • Hard Rock Maple (HM) | • Mahogany (MA)       |                  |
| • Rustic Cherry (RC)   | • Summer Flame (SF)   |                  |

## Weight Rating

- Standard 1000lbs.

## Shipping

The Arezzo Sleeper is shipped assembled.