

Couture Recliner

The Couture Recliner is designed exclusively for healthcare environments, for both treatment applications such as dialysis and for patient rooms.

Standard Features

The Couture Recliner is “fully loaded” with the following standard features:

- Field replaceable back and seat upholstery
- Flow-through design minimizes fluid retention
- Nylon-webbed seat deck and inner back structure provides enhanced comfort for longer procedures
- Heavy duty design throughout for extra strength and product longevity
- Lifetime Warranty on construction; two year warranty on mechanism
- 500lb weight capacity
- Arm cap designed to assist patient ingress and egress
- Infinite recline positioning
- Geometrically balanced mechanism substantially reduces effort required to operate leg extension
- Side lever permits both the patient and caregiver to operate leg extension
- Available in two seat widths: 23” and 26”
- Seat and back foam meets or exceeds CAL 117-2013 flammability requirements



Optional Features

- Swing-up arm for easier patient transfer
- Wood armcap
- Polyurethane armcap
- IV pole w/mounting bracket (can be mounted to either side)
- Central caster locking system
- Trendelenburg positioning
- Flip-up side tray (can be mounted to either side)
- Slide-out foot rest for easier patient transfer
- Moisture Barrier on seat
- CAL TB 133 compliance or certification
- Foley Bag hook
- Upholstered arms
- Wooden legs
- Metal legs

Construction Specifications

- Inner back frame structure is constructed from steel (not plywood). It is lightweight yet very durable
- Heavy-duty long life mechanism tested to 35,000 + cycles
- Arm caps available in 7 attractive finishes shown below:
 - Candlelight (CL)
 - Hard Rock Maple (HM)
 - Summer Flame (SF)
 - Cherry Blossom (CB)
 - Mahogany (MA)
 - Chocolate (CH)
 - Rustic Cherry (RC)