

ergozoneTM

Advanced Pressure Redistribution Mattress



Wound Prevention. Superior Comfort.

ErgozoneTM advanced pressure redistribution mattresses are designed specifically for patient comfort and pressure ulcer prevention. Foam sculpting, dedicated envelopment zones and visco-elastic foam work in unison to maximize envelopment, minimize pressure and increase patient satisfaction.

Ergozone™ – advanced pressure redistribution

Highlights

- Enhanced pressure redistribution zones for high risk areas
- Optimal envelopment
- Top cover helps manage shear and friction
- Increased patient comfort and satisfaction




Pressure Ulcer Prevention

FOAM SCULPTING

creates alterations in the mattress (by pattern, size, depth, spacing and location) to maximize pressure redistribution and protect the most vulnerable areas of the body.

- Greater air flow for temperature and moisture management
- Increased support retention over the life of the mattress



A Portfolio of Prevention & Comfort

	ergozone™ ₃₀₀₀ Performance	ergozone™ ₅₀₀₀ Premium	ergozone™ _{IGEL} Elite
5-Year Foam/Core; 2-Year Cover: Prorated after 1-year	•	•	•
Environmentally Friendly Manufacturing Process with natural ingredients for sustainability	•	•	•
Patient Weight Support up to 500 lbs	•	•	•
Enhanced Pressure Redistribution Zones engineered for high risk areas including the head, shoulders, hips, legs and heels	3 zones	5 zones	infinite zones
Envelopment Zones strategically combine with body contouring zones to enhance pressure redistribution	•	•	•
Visco-Elastic Foam allows greater air flow for temperature and moisture control	•	•	•
4-Way Stretch Cover moves in every direction to minimize shear and friction	•	•	•
Gel Infused Surface with Phase Change Material offers enhanced heat dissipation, moisture management and patient comfort			•