



Overbed Tables - ST140 & ST160 Standards & Specifications

Designed to complement our patient room case good designs, our Overbed Tables combine heavy-duty construction with attractive, functional design.

Tops

Top consists of a 1" thick MDF substrate surfaced with woodgrain thermofoil laminate. Bottom side is finished in white. Foil thickness is .016" and is available in Linet's (7) popular finishes shown below:

- Candlelight (CL)
- Chocolate (CH)
- Hard Rock Maple (HM)
- Mahogany (MA)
- Porto (PO)
- Sugar Maple (SM)
- Summer Flame (SF)

Top dimensions are as follows:

- Model ST140: 16" deep x 30" wide. 3/4" thickness
- Model ST160: 16" deep x 30" wide. 1" thickness

The thermoformed top provides a certain softness, both in appearance and feel. It also eliminates bacteria-collecting edge laminate and features a recessed surface for spill containment. (3/16" deep on Model ST160 and 1/8" deep on Model ST140)

Top is fastened to column by means of (6) wood screws and can be replaced easily in the field.

Base and Column

Base and Column are constructed with 14 gauge steel available in 3 finish choices: black or almond powdercoat; or chrome plated. Column height range is 30" - 44" with 15 locking positions.

Weight Capacity

- Model ST140: 50 lb.
- Model ST160: 100 lb.