



iGel - BodyZone™



Advanced Pressure Redistribution Mattress

Enhanced pressure distribution and optimal envelopment make the iGel mattress superior to the industry standard for prevention support surfaces. Engineered specifically to reduce pressure to high risk areas and to conform to the patient's body, iGel offers excellent patient comfort, which aids in the healing process. Linet's iGel technology maintains an ambient temperature that decreases the patient's need to change positions and prevents tossing and turning.

BodyZone™ – Put your patients in the zone of full body comfort

Prevention Mattress

- VPF™ AND SMT™ TECHNOLOGY GIVE OPTIMAL ENVELOPMENT AND EXCELLENT PRESSURE REDISTRIBUTION PROVIDING SUPPORT FOR PATIENT WEIGHTS UP TO 500LBS
- SMT™ SIGNIFICANTLY REDUCES SHEAR AND FRICTION
- OPEN CELL VISCO TECHNOLOGY PROVIDES ADVANCED TEMPERATURE AND MOISTURE MANAGEMENT WHILE BEING ENTIRELY ANTI-MICROBIAL

- EVERY BODYZONE MATTRESS IS MANUFACTURED USING AN ENVIRONMENTALLY FRIENDLY PROCESS (VPF™) AND USES RENEWABLE MATERIALS
- THE 360° STRETCH COVER* MOVES IN EVERY DIRECTION AND IS IMPERVIOUS TO FLUIDS
- MATTRESSES ARE AVAILABLE IN ALL SIZES AND INCORPORATE PROPRIETARY VPF™ AND SMT™ TECHNOLOGY

VARIABLE PRESSURE FOAMING PROCESS

- VPF™ is a technology used to create different levels of support with high envelopment and maximum comfort. This leads to better height and support retention over the life of the mattress.


SURFACE MODIFICATION TECHNOLOGY

- SMT™ is used to create alterations in pattern, size of pattern, spacing between patterns and depth of cut to protect the most vulnerable areas of the body. This provides a zoned-cushioning layer, enhances air circulation and maximized comfort.

iGEL BODYZONE™ MATTRESS

- Specialized gel technology infused into a honeycomb molded mattress core
- Visco, iGel, and SMT technologies are used in combination to redistribute pressure
- Over 50% less product weight versus other gel based mattresses
- iGel will not leak, dry out, or bottom out

THE LINET ERGOFRAME® IN COMBINATION WITH THE iGEL MATTRESS PROVIDES PREMIUM PATIENT COMFORT BY NEUTRALIZING FRICTION POINTS AND MINIMIZING PRESSURE LOADS. THIS LEADS TO IMPROVED PATIENT OUTCOMES AND QUICKER RECOVERY TIME.



Bodyzone Features & Benefits

A PORTFOLIO OF PRESSURE REDISTRIBUTION MATTRESSES

	BodyZone™ 3000 Premium	BodyZone™ 5000 Premium Plus	BodyZone™ 7000 iGel
360° Cover *	●	●	●
VPF™	●	●	●
SMT™	●	●	●
Natural Ingredients	●	●	●
Open Cell Visco	●	●	●
Microspring™		●	●
True Suspension			●
Infused gel technology			●

Linnet mattresses meet or exceed the requirements of 16 CFR 1632 and 16 CFR 1633.

I - Figures taken from a clinical abstract entitled »Pressure redistribution and envelopment with BodyZone™ mattresses: results of a series of tests.«

II - BODYZONE, SMT and VPF are registered trademarks of Foamex Innovations Operating Company.



10420-R Harris Oaks Blvd, Charlotte, NC 28269 USA, www.linnetamericas.com
Ph: 704-248-5650, Fax: 704-248-5655, e-mail: info@linnetamericas.com