

1st LEVEL of Early Mobilization with LINET

Micro-shifting



Micro-shifting makes it possible to tilt the patient one degree at a time.

Automatic Lateral Therapy



Automatic Lateral Therapy (ALT) is a platform-based Lateral Tilt which can be individually programmed and is recommended for patients unable to tolerate manual turning. ALT can be used in cases of higher risk of VAP which are unable to train the patient's body to tolerate side to side movement.

30° stop



30° stop maintains optimal head of bed (HOB) angle in compliance with pressure injury prevention and Ventilator Associated Pneumonia prevention.

Lateral Tilt



The frame-based Lateral Tilt is a unique LINET feature which can help accelerate recovery of critical care patients.

Discover next level!

Interested in in-depth learning?