

Ave 2 Birthing Bed



Welcome to the world

The journey begins

The innovative design of the AVE 2 brings a new level of comfort and efficiency to every phase of childbirth.

AVE 2 is an ideal platform for safely performing labor exercises with an extra low heigh of 23.6 inches.



LINET | Ave2



Innovation

Simple swivel mechanism allows the foot section to retract and rotate away from the patient quickly and easily.

- No lifting a heavy foot section reduces the chance of caregiver injury.
 Foot section never touches the floor, reducing contamination and trip hazard issues.
 Facilitates rapid and efficient configuration of the bed.
- Enhanced access for caregivers.



The innovative seat sections allows for quick configuration and new positioning options for caregivers and patients.



When tilted, the section can be also used as a foot support.



LINET | Ave2

Simple operation

Rapid configuration with simple and intuitive controls improves efficiency and reduces care team frustration.



Quick adjustment of the leg rests into the required position.



U-SHAPE

The unique oval shape improves perineal access for procedures and minimizes gap around birthing tub.



HAND GRIPS Stable and comfortable grips.

COLUMN DESIGN

Innovative column design allows for both a low height for patient egress and greater height for improved caregiver access during procedures.





Motorized pelvic tilt increases positioning options



Features

Siderails have integrated controls for both patient and caregiver that easily rotate down for patient transfer or access.



SIMPLE CONTROL

Integrated control panel with clear graphics allow for simple operation by patient and caregiver. Lockout of specific functions improves safety.

Controls:

- 1. High and Tilt Adjustment
- 2. Trendelenburg Position
- 3. Lockout Control
- 4. Under Bed Light
- 5. Nurse Call (optional)



SIMPLIFIED CLEANING

High hygienic standard and compact construction of particular sections of the rest area and the vinyl cover of the bed base reduce the time needed for the cleaning of the bed.

- 1. Seamless upholstery.
- 2. Simple mattress removal for easier cleaning access.
- 3. Epoxy powder coating paint for better durability and corrosion resistance.

The smooth design angles, column design and open architecture improve safety by reducing pinch points and sharp edges.







LINET | Ave2

- 1. Trendelenburg
- 2. CPR
- 3. Battery backup
- 4. Easily clean and disinfect
- 5. Lockout control
- 6. Innovative column design
- 7. Easy to access central brake pedal
- 8. Steering caster and handles for better maneuverability
- 9. Removable headboard for anesthesiologist
- 10. Nurse Call (optional)
- 11. Tilting holders for IV pole

Positioning options for

every patient







Semi recumbent position with the feet rested on the foot rests



Semi recumbent position with the feet on the foot section



Side lying





Labor bar



All fours position



Half sitting using foot section



Perineal repair



Half sitting using leg holders



Half sitting with partner assistance

Ave 2 at a glance

COLOR OPTIONS

Color options to complement any room design.



TECHNICAL SPECIFICATIONS

Total length of the lying area	83.54 in	Seat section angle	0°/+18.5°
Minimum length	61.81 in	Trendelenburg angle	0°/-10°
Full width	38.58 in	Leg rests vertical angle	0°/+135°
Castor size	5.9 inch	Leg rests horizontal angle	0°/+60°
Minimum height	23.62 in	Foot section adjustment angle	0°/-22°
Maximum height	41.33 in	Maximum patient load	529 lbs
Back section angle	-10°/+70°	Bed weight	540 lbs



© 2017 LINET Americas, Inc. ALL RIGHTS RESERVED. | LA-0087-REV1