

# Couture Treatment Recliner



## INFINITE RECLINING POSITIONS

Thoughtfully designed and engineered for ease of use in rigorous 24/7 healthcare environments, the Couture Treatment Recliner is also generously padded for all over patient comfort.

## FEATURES

- Field replaceable back and seat upholstery
- Flow-through design minimizes fluid retention
- Nylon-webbed seat deck and inner back structure provides enhanced comfort for longer procedures
- Heavy duty design throughout for extra strength and product longevity
- Arm cap designed to assist patient ingress and egress
- Weight Rating: Standard 500 lb
- Geometrically balanced mechanism reduces effort required to operate leg extension
- Side lever permits both the patient and caregiver to operate leg extension
- Available in two seat widths: 23" and 26"
- Seat and back foam meets or exceeds CAL 117-2013 flammability requirements
- One year parts warranty

# COUTURE OPTIONS



*Shown with optional wood arm cap, central locking system and Trendelenburg feature. Cover image also features optional IV pole, side tray and foley bag holder.*


## Additional Options

- Swing-up arm for easier patient transfer
- IV pole with mounting bracket
- Central castor locking system
- Trendelenburg positioning
- Polyurethane armcaps in black, grey or taupe
- Flip-up side tray
- Slide-out foot rest for easier patient transfer
- Moisture barrier on seat
- CAL TB 133 compliant or certified
- Foley bag hook

Wood armcaps are available in these standard wood finishes:



Hardrock  
Maple



Candlelight



Rustic Cherry



Summer Flame



Mahogany



Cherry  
Blossom



Grenada



Chocolate

Visit [www.LINETAmericas.com](http://www.LINETAmericas.com) for full list of fabric options.



Members of LINET Group

## LINET Americas, Inc.

Address: 10420-R Harris Oaks Blvd., Charlotte, NC 28269

Phone: 877-815-9895 | Fax: 704-248-5655 | E-mail: [info@linetamericas.com](mailto:info@linetamericas.com) | [www.linetamericas.com](http://www.linetamericas.com)

© 2020 LINET Americas, Inc. | 040820  
LINET® is a registered trademark of LINET Group SE