

Ergozone®



WOUND PREVENTION. SUPERIOR COMFORT.

Ergozone® 3000 and 5000 are advanced pressure redistribution mattresses designed specifically for patient comfort and pressure injury prevention. Foam sculpting, dedicated envelopment zones and viscoelastic foam work in unison to maximize envelopment, minimize pressure and increase patient satisfaction.

FOAM SCULPTING

Foam sculpting creates alterations in the top layer of the mattress (by pattern, size, depth, spacing, and location) to maximize pressure redistribution and protect the most vulnerable areas of the body.

- Greater air flow for temperature and moisture management.
- Increases support retention over the life of the mattress.



- Enhanced pressure redistribution zones for high risk areas
- Optimal envelopment for patient comfort
- 5° heel slope
- Top cover helps manage shear and friction
- Increased patient comfort and satisfaction

Ergozone[®] at a glance

A Portfolio of Prevention and Comfort

	Ergozone 3000 Performance	Ergozone 5000 Premium
Environmentally Friendly Manufacturing Process with natural ingredients for sustainability.	●	●
Patient Weight Support up to 500 pounds.	●	●
Enhanced Pressure Redistribution Zones engineered for high risk areas including the head, shoulders, hips, legs, and heels.	3 zones	5 zones
Envelopment Zones strategically combine with body contouring zones to enhance pressure redistribution.	●	●
Open Cell Viscoelastic Foam to allow greater air flow for temperature and moisture control.	●	●
4-Way Stretch Cover moves in every direction to minimize shear and friction.	●	●

TECHNICAL PARAMETERS*

Width	34.25 - 39.5 in	Weight Capacity	500 lb
Length	78.5 -84 in	Heel Slope	5°
Height	5.5 - 6 in	Flammability	Meets Cal TB 129, 16 CFR 1633

**Width, length and height may depend on bed frame model. Compatible LINET[®] bed frames include Eleganza 3, Eleganza 5, Multicare[®], Multicare LE, and Sentida.*



Members of LINET Group

LINET Americas, Inc.

Address: 10420-R Harris Oaks Blvd., Charlotte, NC 28269

Phone: 877-815-9895 | Fax: 704-248-5655 | E-mail: info@linetamericas.com | www.linetamericas.com

© 2020 LINET Americas, Inc. | 040720

LINET[®], Multicare[®], and Ergozone[®] are registered trademarks of LINET Group SE