



# Multicare<sup>®</sup> LE



Designed as an acuity adaptable bed with innovative features for caregiver and patient safety.

Use the innovative standard features to meet the needs of a variety of care settings.



## Lateral Turn Assist

Utilize the unique frame-based lateral turn to assist with routine care to minimize caregiver injury.

## Early Mobilization

Incorporate standard features for use in Early Mobilization Programs and physical therapy.



# Multicare<sup>®</sup> LE at a glance

## STANDARD FEATURES

### Lateral Turn Assist

Unique frame-based turn up to 15°

### Mobi-Lift<sup>®</sup>

Built-in sit-to-stand device for patient to raise and lower the bed

### i-Brake<sup>®</sup>

Intelligent brakes automatically engage after 60 seconds when bed is plugged into power

### Integrated LED Light

Bright LED light in patient handset

### i-Drive

Automatic 5th wheel enhances maneuverability during transport

### Ergoframe<sup>®</sup>

Reduces shear and friction forces caused by patient migration to the foot of the bed

### i-Drive Power<sup>®</sup> (optional)

Motorized patient transport

### One Touch Controls

- 30° HOB
- Cardiac Chair
- CPR
- Autocontour
- Examination position
- Bed egress position
- Trendelenburg
- Vascular knee adjustment
- Electric bed extension

### Hands-Free Foot Controls

Available on both sides of the bed:

- Bed height
- Exam position
- Left and right turn assist



## TECHNICAL PARAMETERS

Safe working load	550 lb	Frame weight	496 lb
Frame height - low*	17.3 in	Trendelenburg angle	13°
Frame height - high	32.4 in	Reverse Trendelenburg	-16°
Frame width	41.3 in	Max knee angle	35°
Frame length extended	93.3 in	Max HOB angle	70°
Frame length retracted	84.6 in		

\*A lower egress height may be obtained by utilizing the turn function.



Members of LINET Group

### LINET Americas, Inc.

Address: 10420-R Harris Oaks Blvd., Charlotte, NC 28269

Phone: 877-815-9895 | Fax: 704-248-5655 | E-mail: [info@linetamericas.com](mailto:info@linetamericas.com) | [www.linetamericas.com](http://www.linetamericas.com)

© 2020 LINET Americas, Inc. | 040720

LINET<sup>®</sup>, Multicare<sup>®</sup>, i-Brake<sup>®</sup>, i-Drive Power<sup>®</sup>, Ergoframe<sup>®</sup>, and Mobi-Lift<sup>®</sup> are registered trademarks of LINET Group SE