

# Overbed Tables

LINET® offers two lines of overbed tables to meet your facility's needs. Create a customized overbed table with the modern Avalon series or choose from the economic Solido line for a functional solution. LINET overbed tables are built with an ultra low base for easy mobility.

## AVALON



The Avalon overbed table offers a range of design options that can be customized to meet your needs, including thermofoil tops with a variety of shapes and finishes. Accessories are available for customization and to further enhance infection control initiatives.

See reverse side for Avalon customization options.

## SOLIDO



The economic Solido overbed table is designed for modern healthcare environments with a laminate top and durable metal base. Customize for patient needs with a fixed or tiltable eating board.

The Solido overbed table does not offer additional customizations.

# Avalon Options

## TABLE TOP SHAPES

Standard: 19" x 32"  
Extended Length: 18" x 36"



Kidney  
(Thermofoil)




Peanut  
(Thermofoil)



Jelly Bean  
(Thermofoil)

Standard: 17" x 32"  
Extended Length: 18" x 36"



Laminate or  
Thermofoil

Standard: 17" x 32"  
Extended Length: Custom



White Board  
(Laminate)

## TABLE TOP FINISHES



Hard Rock  
Maple



Candlelight



Amber Cherry



Acajou  
Mahogany



Finesse  
(Thermofoil only)



Chocolate Pear

## BASE FINISHES



Graphite Gray



Almond



Nevada

## ADDITIONAL ACCESSORIES



Integrated  
Accessory Rail



Low Profile Base Options: U or H shape  
(1.75" clearance)



Ultra Low Profile  
Base (1.5" clearance)



Storage Baskets  
9x9x10, 10x6x9, 5x5x4



Integrated Cup  
Holders



Flip Surface



Hand Sanitizer



Glove Box/  
Tissue Holder



Dual Vanity Drawer  
and Mirror



Chart Holder



Urinal  
Hook



Cord  
Wrap



Power  
Bar



Members of LINET Group

### LINET Americas, Inc.

Address: 10420-R Harris Oaks Blvd., Charlotte, NC 28269

Phone: 877-815-9895 | Fax: 704-248-5655 | E-mail: [info@linetamericas.com](mailto:info@linetamericas.com) | [www.linetamericas.com](http://www.linetamericas.com)

© 2020 LINET Americas, Inc. | 040820

LINET® is a registered trademark of LINET Group SE